# Física Experimental II

## uniandes

Departamento de Física

2015

Guía de laboratorio

**Física Experimental II**

*Guía de laboratorio*

Revisión 2015:

Juan David Lizarazo Ferro

Benjamín Oostra

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# Manual del calorímetro

## 1. PURPOSE AND DESCRIPTION

Calorimeter 04401.00 is a water calorimeter with which the specific heat of solids or liquids can be determined. Conversion energies, such as the melting energy of ice, can also be determined.

The unit is fitted with an electric heating system to heat up the contents. Power is supplied to the heating resistor over two 4 mm sockets in the lid. The lid has an orifice ( $d = 10$  mm) to introduce a thermometer or a temperature probe. An agitator with a yoke handle, which can be lifted, passes through two smaller orifices in the lid. The calorimeter is designed to stand temperatures up to 100 °C.

## 2. HANDLING

The handling is explained taking the determination of specific thermal capacities as an example. In the case of solids, the mixing method is used, for liquids, energy is added by means of the electric heating element.

In order to accelerate thermal equilibrium in the calorimeter, the agitator should be continuously moved during the experiment. It may, however, not be lifted to the point where water is splashed onto the styrofoam lining of the lid.

### 2.1 Mixing metod

To determine the specific heat of solids according to the mixing method, a sample body of known temperature and mass is brought into thermal contact with a quantity of water of known temperature and thermal capacity inside the calorimeter. The specific heat of the sample is calculated from the temperature of the mixture after thermal equilibrium is reached.

#### *Carrying out of experiment:*

- An adequate quantity of liquid (200 ml or more) is weighed and filled into the calorimeter.
- The sample body is hanged from a support by means of a fishing thread and heated in a sufficiently stationary water bath, e. g. in boiling water to 100 °C. A gauze bag may be used to carry several small samples.
- The temperature of the heated sample (temperature of the bath) and the temperature of the water in the calorimeter are read as precisely as possible from the corresponding thermometers immediately before immersing the sample in the calorimeter.
- When one is certain the sample body has reached the temperature of the bath, it is immersed as fast as possible in the calorimeter. The calorimeter lid is closed at once and agitation starts.
- When temperature in the calorimeter begins to drop due to thermal release to the environment, maximum temperature is read as temperature of the mixture.

### 2.2 Electric heating method

To determine the specific heat of liquids, a quantity of liquid of known mass and temperature is filled into the calorimeter and heated by the electric heating element. The specific heat of the liquid is calculated from energy input and temperature increase.

Next to water, only such liquids may be filled into the calorimeter, which do not attack aluminium, nickel or styrofoam, e. g. all types of alcohol. Only alternating current should be used, to avoid corrosion of the heating element. It is furthermore recommended to use highly purified (distilled) water.

The heating element can only be operated when it is completely immersed in the liquid. A filling quantity of 200 ml is sufficient for this.

The amount of added electric energy is determined measuring current intensity, voltage and heating time. Adequate current intensities: 3 ... 5 A. Adjustment of the supply voltage should be determined in a preliminary experiment, so that the required heating energy will be immediately available during the main experiment.

*Carrying out the experiment:*

- An adequate amount of liquid (200 ml or more) are weighed and filled into the calorimeter.
- Initial temperature is read; voltage supply and chronometer are switched on; continuous agitation is assured.
- After temperature has increased by 5 - 10 °C for example, voltage supply and chronometer are switched off simultaneously and the maximum value displayed by the thermometer, before temperature stops increasing, is read.

To obtain a more precise measurement, thermal energy released by the calorimeter to the environment must be taken into account. This is achieved by means of a correction of the read final temperature:

- a second experiment is carried out with the same calorimeter contents, during which the drop of temperature  $\epsilon$  is measured at the average heating up temperature during a period of time which corresponds to the total heating time during the main experiment. Calculation is now repeated with the final temperature increased by  $\epsilon$ .

### 3. EXPERIMENTING REFERENCE LITERATURE

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OPERATING MANUAL

## **MODEL YF-3503**

I. Thanks for your patronage. Read this operating manual before beginning. Correct operations help develop the best performances of multimeter also decrease the possibility of damages.

### **II. Specification**

#### **2-1 General Specification**

- ( 1) LCD: 3½ digits, with a maximum reading of 1999.
- ( 2) Polarity Indication: Automatic polarity, " - " display for negative input.
- ( 3) Overload Indication: LCD will show a "1" or " - 1" in the left highest position.
- ( 4) Low Battery Indication: Replace battery as LCD display "BAT".
- ( 5) Sample-Rate: 2.5 times per second.
- ( 6) Power Supply: 006P DC9V, 1 pc.
- ( 7) Battery Life: 200 hours approx.
- ( 8) Operating Temperature & Humidity: 0°C – 40°C below 80% RH.
- ( 9) Dimension & Weight: 143 × 74 × 38 mm, 257g approx. (with Battery)
- (10) Accessory: Operating manual, test leads, battery 006P DC9V 1pc.

#### **2-2 Electrical Specifications (23°C ± 5°C, below 80% RH) Accuracy: '± (..... % rdg + ..... dgt)**

#### **DCV -**

Range	Resolu-tion	Accuracy	Input Impedance	Overload Protection
200mV	0.1mV	0.5% + 1		DC 500V AC 350V rms
2V	1mV			
20V	10mV			
200V	100mV	0.8% + 1	10MΩ	DC 1100V AC 800V rms
1000V	1V			

### ACV

Range	Resolu-tion	Accuracy	Input Impedance	Overload Protection
200mV	0.1mV	1.2% + 3 (40 Hz - 500 Hz)	10MΩ	DC 500V AC 350V rms
2V	1mV			DC 1100V AC 800V rms
20V	10mV			
200V	100mV			
750V	1V			

### DCA

Range	Resolu-tion	Accuracy	Burden Voltage	Overload Protection
200μA	0.1μA	0.8% + 1	0.3V	0.5A Fuse & Diodes
2mA	1μA		0.3V	
20mA	10μA		0.3V	
200mA	100μA		0.3V	
20A	10mA	2% + 1	0.7V	

### ACA

Range	Resolu-tion	Accuracy	Burden Voltage	Overload Protection
200μA	0.1μA	1.2% + 3 (40 Hz - 500 Hz)	0.3V max	0.5A Fuse & Diodes
2mA	1μA		0.3V max	
20mA	10μA		0.3V max	
200mA	100μA		0.3V max	
20A	10mA	2% + 3	0.7V max	

### Resistance ( Ω )

Range	Resolu-tion	Accuracy	Open Voltage	Overload Protection
200Ω	0.1Ω	1% + 2	± 2.8V	DC/AC 500V rms
2KΩ	1Ω	0.8% + 2	± 0.35V	
20KΩ	10Ω		0.35V	
200KΩ	100Ω		0.35V	
2MΩ	1KΩ		0.35V	
20MΩ	10KΩ	2% + 2	0.35V	

## Diode Test

Range	Resolution	Test Current	Open Voltage	Overload Protection
→	1mV	1mA	≤ 2.8V	DC/AC 500 V rms

## Continuity Test •11)

Range	Resolution	Sound	Open Voltage	Overload Protection
•11)	1 Ω	below 40Ω	below ≤ 2.8V	AC/DC 500 V rms

## hFE

Range	Range	Test Voltage	Test Current	Overload Protection
NPN PNP	0-1000	< 3.2V	IB ≤ 10μA	

## Battery

Range	Resolution	Load Current	Load Resistance
1.5V	1mV	80mA	19Ω
9V	10mV	7mA	1.3KΩ

## LOGIC TTL


Range	Input Impedance	Input Voltage	Logic Level	Min detect impulse	Max output frequency	Overload protection
TTL ↓	120K Ω	DC5V	I 2.8 ± 0.8V *▲* 0 0.8 ± 0.5V *▼*	25nS	20MHz	500Vrms

## CAPACITANCE

Range	Resolution	Accuracy	TEST FREQUENCY
2nF	1PF		
20nF	10PF		
200nF	100PF	3% + 10	
2μF	1aF		400 Hz
2μF	1nF		

**III. Name of parts and operation (as the picture below)**

- (1) LCD: 3½ digits with a maximum reading of 1999, also including indications of " - ", low battery.
- (2) AC/DC Slide Switch: For selection of "AC" or "DC" measurements.
- (3) Range Rotary Selector: For power on and off, and selection of desired function range.
- (4) 20A Measuring Connector: For measuring a current below 20A.
- (5) COM Measuring Connector: To connect negative lead (black test lead) for all kinds of measurements.
- (6) VΩA Measuring Connector: To connect positive lead (red test lead) for voltage, current, resistance measurements.
- (7) hFE Measuring Connector: To insert transistor for hFE measurement.
- (8) Capacity Measuring Connector.
- (9) Power Switch: ON-Power ON.  
OFF-Power OFF.

#### **IV. Measurement**

##### **4-1. DCV (ACV) Measurement**

- (1) Plug red test lead to "VΩA" connector, and black test lead to "COM" connector.
- (2) Place switch to "DC" or "AC" function according to the need.
- (3) Turn the function range selector to "V" region (if measuring an unknown voltage, start from the highest range, then adjust to a proper lower range for a best resolution).
- (4) Connect the other ends of test leads to the desired circuit.
- (5) Read the value from the LCD.

##### **4-2. DCA (ACA) Measurement**

- (1) Plug red test lead to "VΩA" connector (if measuring a current between 200mA and 20A, plug to "20A" connector instead), and black test lead to "COM" connector.
- (2) Place switch to "DC" or "AC" function according to the need.
- (3) Turn the function range selector to "A" region (if measuring an unknown current, start from the highest range, then adjust to a proper lower range for a best resolution).
- (4) Connect the other ends of test leads in series to the desired circuit.
- (5) Read the value from the LCD.

##### **4-3. Resistance Measurement**

- (1) Plug red test lead to "VΩA" connector, and black test lead to "COM" connector.
- (2) Turn the function range selector to "Ω" region (if measuring an unknown resistance, start from the highest range, then adjust to a proper lower range for a best resolution).
- (3) Connect the other ends of test leads to the desired resistor.
- (4) Read the value from the LCD.

##### **4-4. Diode $\left(-\right)$ Measurement**

- (1) Turn the function range selector to " $\left(-\right)$ " region.
- (2) Plug red test lead to "VΩA" connector, and black test lead to "COM" connector.

- (3) Connect the other end of red test leads to the positive pole (P) of the diode, and the other end of black test leads to the negative pole (N).
- (4) Read the value from the LCD. (If connection to poles of diode is reverse, LCD will display "1" in the highest digit.)

#### 4-5. Continuity (-■) Measurement

- (1) Turn the function range selector to "-■" region.
- (2) Plug red test lead to "VΩA" connector, and black test lead to "COM" connector.
- (3) Connect the other ends of test leads to test the conduction situation of circuit. If impedance in circuit is below 40Ω, the buzzer will sound.

#### 4-6. hFE Measurement (for Transistor)

- (1) Turn the function range selector to "hFE" (NPN or PNP) region.
- (2) Insert the transistor pins into "hFE" connector according to E.B.C. sequence.
- (3) Read the value from LCD.

#### 4-7. Battery Power (BAT.) Measurement

- (1) Plug red test lead to "VΩA" connector, and black test lead to "COM" connector.
- (2) Turn the function range selector to "1.5V" or "9V" range.
- (3) Connect the other end of test leads to the desired battery, red test lead to positive pole, and black test lead to negative pole.
- (4) Read the value from the LCD.

#### 4-8. Logic Measurement

- (1) Turn the function switch on the "TTL" position.
- (2) Connect the red test lead to the "VΩA" connector, the black to the "COM" connector.
- (3) Confirming the logic circuit voltage at 5 V level, connect the black test lead to the logic circuit the negative end of the electric source, with the red test lead do logic circuit test.
- (4) When the logic level is "1", LCD will display the "▲" mark. When the logic level is "0", LCD will display "▼" mark, and the buzzer will sound.

#### **4-8. Capacity (C) Measurement**

- (1) Turn the function range selector to "F" region (if measuring an unknown capacitor, start from the highest ranges, then adjust to a proper lower range for a best resolution).
- (2) Insert the desired capacitor into the measuring connector.
- (3) Read the value from the LCD.

### **V. Maintenance**

#### **5-1. Battery Replacement**

As battery power is not sufficient, LCD will display "BAT", replace battery to keep measurement accurate.

- (1) Remove all the test leads, and turn the power off.
- (2) Open the battery cover according to instructed direction, take out the battery, and install a new one.
- (3) Secure the battery cover.

#### **5-2. Fuse Replacement**

If replacement is required, please replace with a same one rated 0.5A 250VAC, 5 $\times$  20mm to keep circuit protected normally.

- (1) Turn the power off.
- (2) Open the battery cover according to instructed direction, take out the burned fuse, and install a new one.
- (3) Secure the battery cover.

#### **5-3. Note**

This multimeter is a precision instrument, do not violate the specifications no matter in use or in storage. Besides, following circumstances should be avoided to keep measurement accurate.

- (1) Do not measure or connect DCV over 1000V, or ACV over 750V.
- (2) As function range selector sets on " $\Omega$ ", "hFE", " $\mu A$ ", ".1") " ranges, do not abb any source of voltage or current.
- (3) Take out battery for long time storage, also avoid storing in places in high temperature and humidity.

# User Manual

## APÉNDICE C

---

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---

**PS280 & PS283  
DC Power Supplies**

**070-8355-03**

Copyright © Tektronix, Inc. 1991. All rights reserved.

Tektronix products are covered by U.S. and foreign patents, issued and pending. Information in this publication supercedes that in all previously published material. Specifications and price change privileges reserved.

Tektronix, Inc., P.O. Box 1000, Wilsonville, OR 97070-1000

TEKTRONIX and TEK are registered trademarks of Tektronix, Inc.

## **WARRANTY**

Tektronix warrants that this product will be free from defects in materials and workmanship for a period of one (1) year from the date of shipment. If any such product proves defective during this warranty period, Tektronix, at its option, either will repair the defective product without charge for parts and labor, or will provide a replacement in exchange for the defective product.

In order to obtain service under this warranty, Customer must notify Tektronix of the defect before the expiration of the warranty period and make suitable arrangements for the performance of service. Customer shall be responsible for packaging and shipping the defective product to the service center designated by Tektronix, with shipping charges prepaid. Tektronix shall pay for the return of the product to Customer if the shipment is to a location within the country in which the Tektronix service center is located. Customer shall be responsible for paying all shipping charges, duties, taxes, and any other charges for products returned to any other locations.

This warranty shall not apply to any defect, failure or damage caused by improper use or improper or inadequate maintenance and care. Tektronix shall not be obligated to furnish service under this warranty a) to repair damage resulting from attempts by personnel other than Tektronix representatives to install, repair or service the product; b) to repair damage resulting from improper use or connection to incompatible equipment; or c) to service a product that has been modified or integrated with other products when the effect of such modification or integration increases the time or difficulty of servicing the product.

**THIS WARRANTY IS GIVEN BY TEKTRONIX WITH RESPECT TO THIS PRODUCT IN LIEU OF ANY OTHER WARRANTIES, EXPRESSED OR IMPLIED. TEKTRONIX AND ITS VENDORS DISCLAIM ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. TEKTRONIX' RESPONSIBILITY TO REPAIR OR REPLACE DEFECTIVE PRODUCTS IS THE SOLE AND EXCLUSIVE REMEDY PROVIDED TO THE CUSTOMER FOR BREACH OF THIS WARRANTY. TEKTRONIX AND ITS VENDORS WILL NOT BE LIABLE FOR ANY INDIRECT, SPECIAL, INCIDENTAL, OR CONSEQUENTIAL DAMAGES IRRESPECTIVE OF WHETHER TEKTRONIX OR THE VENDOR HAS ADVANCE NOTICE OF THE POSSIBILITY OF SUCH DAMAGES.**

# Table of Contents

<b>General Safety Summary .....</b>	<b>iii</b>
<b>Getting Started .....</b>	<b>1</b>
Preparing the Power Supply for Use .....	2
Front Panel .....	4
Turning On the Instrument .....	8
<b>Reference .....</b>	<b>9</b>
Constant Voltage/Constant Current Crossover .....	9
Setting the Current Limit .....	10
Test Modes .....	11
Independent Modes .....	11
Tracking Modes .....	20
<b>Appendix A: Specifications .....</b>	<b>25</b>
<b>Appendix B: Maintenance .....</b>	<b>29</b>
Cleaning .....	29
Preparing for Shipment .....	29
Troubleshooting .....	30
<b>Appendix C: Replaceable Parts .....</b>	<b>33</b>
Standard Accessories .....	33
Optional Accessories .....	33

# List of Figures

Figure 1: Line Voltage Selectors, Power Input, and Fuse Locations .....	2
Figure 2: PS280 or PS283 Front Panel .....	4
Figure 3: Constant Voltage/Constant Current Crossover .....	10
Figure 4: Independent Floating Application .....	12
Figure 5: Independent Common Ground-Referenced Application .....	14
Figure 6: Independent Ground-Referenced Split Application ..	15
Figure 7: Three Ground-Referenced Negative Power Supplies ..	16
Figure 8: Three Ground-Referenced Positive Power Supplies ..	17
Figure 9: Independent Positive Stacked Application .....	18
Figure 10: Independent Negative Stacked Application .....	19
Figure 11: Series Tracking Inside the PS280 or PS283 .....	20
Figure 12: Series Tracking Application .....	21
Figure 13: Parallel Tracking Inside the PS280 or PS283 .....	23
Figure 14: Parallel Tracking Application .....	24

# General Safety Summary

Review the following safety precautions to avoid injury and prevent damage to this product or any products connected to it.

## Injury Precautions

### **Use Proper Power Cord**

To avoid fire hazard, use only the power cord specified for this product.

### **Avoid Electric Overload**

To avoid electric shock or fire hazard, do not apply a voltage to a terminal that is outside the range specified for that terminal.

### **Ground the Product**

This product is grounded through the grounding conductor of the power cord. To avoid electric shock, the grounding conductor must be connected to earth ground. Before making connections to the input or output terminals of the product, ensure that the product is properly grounded.

### **Do Not Operate Without Covers**

To avoid electric shock or fire hazard, do not operate this product with covers or panels removed.

### **Use Proper Fuse**

To avoid fire hazard, use only the fuse type and rating specified for this product.

### **Do Not Operate in Wet/Damp Conditions**

To avoid electric shock, do not operate this product in wet or damp conditions.

### **Do Not Operate in Explosive Atmosphere**

To avoid injury or fire hazard, do not operate this product in an explosive atmosphere.

## **Product Damage Precautions**

### **Use Proper Voltage Setting**

Before applying power, ensure that the line selector is in the proper position for the power source being used.

### **Provide Proper Ventilation**

To prevent product overheating, provide proper ventilation.

### **Do Not Operate With Suspected Failures**

If you suspect there is damage to this product, have it inspected by qualified service personnel.

## **Safety Terms and Symbols**

### **Terms in This Manual**

These terms may appear in this manual:

---

***WARNING.*** *Warning statements identify conditions or practices that could result in injury or loss of life.*

---

**CAUTION.** Caution statements identify conditions or practices that could result in damage to this product or other property.

### Terms on the Product

These terms may appear on the product:

DANGER indicates an injury hazard immediately accessible as you read the marking.

WARNING indicates an injury hazard not immediately accessible as you read the marking.

CAUTION indicates a hazard to property including the product.

### Symbols on the Product

The following symbols may appear on the product:

DANGER  
High Voltage

Protective Ground  
(Earth) Terminal

ATTENTION  
Refer to  
Manual

Double  
Insulated

## Certifications and Compliances

### CSA Certified Power Cords

CSA Certification includes the products and power cords appropriate for use in the North America power network. All other power cords supplied are approved for the country of use.

# Getting Started

The Tektronix PS280 or PS283 Laboratory DC Power Supply is a multifunction bench or portable instrument. This regulated power supply provides a fixed 5 V output for powering logic circuits and two variable outputs for a wide variety of test and experimental uses.


The PS280 or PS283 contains two identical, independently adjustable DC power supplies that you can vary from 0 to 30 V. The current on the PS280 variable power supplies varies from 0 to 2 A. The current on the PS283 variable power supplies varies from 0 to 1 A. In all other respects the instruments are identical. Unless otherwise noted, descriptions and procedures in this manual apply to both instruments.

Front panel switches select one of three modes of operation:

- Independent — In this mode, the output voltage and current of each supply can be controlled independently.
- Series — In this tracking mode, the variable outputs are connected in series, and the controls of the master power supply adjust the voltages or currents of both power supplies. Series mode allows the power supplies to be varied from 0 to 60 V at 0 to 2 A for the PS280, or 0 to 1 A for the PS283.
- Parallel — In this tracking mode, the variable outputs are connected in parallel, and the controls of the master power supply adjust the voltages or currents of both power supplies. Parallel mode allows the power supplies to be varied from 0 to 4 A for the PS280 at 0 to 30 V, or from 0 to 2 A at 0 to 30 V for the PS283.

## Preparing the Power Supply for Use

Check the following items prior to operating the Laboratory DC Power Supply for the first time (see Figure 1 for locations of items 1 through 3):

**Figure 1: Line Voltage Selectors, Power Input, and Fuse Locations**

**CAUTION.** To prevent damage to the instrument, set the line voltage selectors to the proper voltage setting and install the correct line voltage fuse before operating the equipment.

1. Set the line voltage selectors to the input line voltage. These selectors connect internal wiring for various line voltages. This product is intended to operate from a power source that does not supply more than 250 V<sub>RMS</sub> between the supply conductors or between either supply conductor and ground. For line voltage ranges, refer to *Appendix A: Specifications* on page 25.

**WARNING.** To prevent electrical shock, unplug the power cord and disconnect the test leads from the circuit before checking or replacing the fuse.

2. Check that the correct line fuse is installed. The line fuse provides protection if the equipment malfunctions or an overload

occurs. Refer to *Appendix C: Replaceable Parts* on page 33 for fuse part numbers.

***WARNING.*** *To prevent electrical shock, connect the power cord to a properly grounded power source. The outside (ground) of this connector is connected through the equipment to the power source ground. Do not remove the ground lug from the power cord for any reason.*

---

3. Connect the input power cord. Use only the power cords specified for this equipment. Refer to *Appendix C: Replaceable Parts* on page 33 for power cord part numbers.

## Front Panel

Figure 2 shows the front-panel controls, connectors, and indicators with brief descriptions following the figure.


Figure 2: PS280 or PS283 Front Panel

1. LED Display. Lights when the instrument is turned on. The numbers indicate the voltage or current produced by the left variable power supply.
2. AMPS/VOLTS Switch. This switch selects whether the LED display for the left variable power supply shows the current or the voltage. If the switch is pushed to the left, the display shows the current. If the switch is pushed to the right, the display shows the voltage.

3. AMPS Indicator. Lights when AMPS is selected with the AMPS/VOLTS switch for the left variable power supply.
4. VOLTS Indicator. Lights when VOLTS is selected with AMPS/VOLTS switch for the left variable power supply.
5. AMPS Indicator. Lights when AMPS is selected with the AMPS/VOLTS switch for the right variable power supply.
6. VOLTS Indicator. Lights when VOLTS is selected with AMPS/VOLTS switch for the right variable power supply.
7. AMPS/VOLTS Switch. This switch selects whether the LED display for the right variable power supply shows the current or the voltage. If the switch is pushed to the left, the display shows the current. If the switch is pushed to the right, the display shows the voltage.
8. LED Display. Lights when the instrument is turned on. The numbers indicate the voltage or current produced by the right variable power supply.
9. POWER Button. Turns on the instrument when pressed. When pressed again, it turns off the instrument.
10. CURRENT Knob. Use this control to set the output current for the right, variable power supply. If the instrument is in a tracking mode, the left power supply is the slave and the CURRENT knob has no effect.
11. C.C. Indicator. If this is lighted, the left variable power supply is producing a constant current. See Figure 3 on page 10 for an illustration of the constant voltage/constant current cross-over point.
12. C.V. Indicator. If this is lighted, the left variable power supply is producing a constant voltage. See Figure 3 on page 10 for an illustration of the constant voltage/constant current cross-over point.

**13. Output Terminals.** These terminals for the left, variable power supply allow you to plug in the test leads as follows:

- The red terminal on the right is the positive polarity output terminal. It is indicated by a + sign above it.
- The black terminal on the left is the negative polarity output terminal. It is indicated by a – sign above it.
- The green terminal in the middle is the earth and chassis ground.

**14. VOLTAGE Knob.** Allows you to set the output voltage for the left variable power supply. If the instrument is in a tracking mode, the left power supply is the slave and the VOLTAGE knob has no effect.

**15. TRACKING Buttons.** These buttons select the test mode of the instrument. The PS280 or PS283 features two tracking modes: series and parallel. If both push-button switches are disengaged (out), the two variable power supplies operate independently. If the left switch is pushed in, the instrument operates in series mode. If both switches are pushed in, the instrument operates in parallel mode.

In series mode, the master power supply controls the voltage for both power supplies, which can then range from 0 to 60 V. Refer to *Series* on page 20 for further details.

In parallel mode, the master power supply controls both the voltage and the current for both power supplies. The current can then range from 0 to 4 A (0 to 2 A for the PS283). Refer to *Parallel* on page 22 for further details.

**16. CURRENT Knob.** Use this control to set the output current for the right, variable power supply. If the instrument is in a tracking mode, the right power supply is the master and the CURRENT knob affects both variable power supplies.

**17. Output Terminals.** These terminals for the right, variable power supply allow you to plug in the test leads as follows:

- The red terminal on the right is the positive polarity output terminal. It is indicated by a plus (+) sign above it.
- The black terminal on the left is the negative polarity output terminal. It is indicated by a minus (-) sign above it.
- The green terminal in the middle is the earth and chassis ground.

**18. C.C. Indicator.** If this is lighted, the power supply is producing a constant current. See Figure 3 on page 10 for an illustration of the constant voltage/constant current crossover point.

**19. C.V. Indicator.** If this is lighted, the power supply is producing a constant voltage. See Figure 3 on page 10 for an illustration of the constant voltage/constant current crossover point.

**20. VOLTAGE Knob.** Allows you to set the output voltage for the right variable power supply. If the instrument is in a tracking mode, the right power supply is the master and the VOLTAGE knob affects both variable power supplies.

**21. Output Terminals.** These terminals for the 5 V FIXED power supply allow you to plug in the test leads as follows:

- The red terminal on the right is the positive polarity output terminal.
- The black terminal on the left is the negative polarity output terminal.

**22. The overload indicator** lights when the current on the 5 V FIXED power supply becomes too large.

## Turning On the Instrument

After you have ensured that the PS280 or PS283 is set up for the proper line voltage and has the proper fuse (refer to *Preparing the Power Supply for Use* on page 2), you are ready to turn it on.

**CAUTION.** *To avoid damaging the PS280 or PS283, do not use it when the ambient air temperature exceeds 40° C. Also, allow adequate space at the rear of the instrument to permit the heat sink to radiate heat.*

---

1. Ensure that the **POWER** button is disengaged (out) and that the instrument is turned off.
2. Plug the power cord into an appropriate power source.
3. Turn both **VOLTAGE** knobs counterclockwise to the minimum setting.
4. Press the **POWER** button. The LED displays light up.

# Reference

This section tells how to set the PS280 or PS283 current limit. It also explains the constant voltage/constant current crossover characteristic of the instrument. Finally, the section includes procedures for using the instrument in both independent and tracking modes and provides examples of a variety of applications.

## Constant Voltage/Constant Current Crossover

The PS280 or PS283 DC Power Supply features a constant voltage/constant current automatic crossover. This feature permits continuous operation in the transition from constant-voltage mode to constant-current mode as the load changes. The intersection of the constant-current and constant-voltage modes is called the crossover point. Figure 3 on page 10 shows the relationship between the load and the crossover point.

For example, if the load is such that the power supply is operating in constant-voltage mode, the PS280 or PS283 provides a regulated output voltage. The output voltage remains constant as the load increases until the preset current limit is reached. Then the crossover occurs. At that point, the output current becomes constant and the output voltage drops in proportion to further load increases.

Crossover is indicated by the front panel red C.C. and green C.V. indicator lights. If the C.V. indicator is lighted, the instrument is operating in constant-voltage mode. If the C.C. indicator is lighted, the instrument is operating in constant-current mode.

Crossover from the constant-current mode to the constant-voltage mode also occurs automatically in response to a decrease in load. For example, suppose you are charging a 12 V battery. Initially, the open circuit voltage of the power supply is preset for 13.8 V. A low battery places a heavy load on the power supply, and it operates in constant-current mode. You adjust the instrument to charge the battery at the rate of 1 A. As the battery becomes charged and its voltage approaches 13.8 V, the load decreases to the point where the battery no longer demands the full 1 A charging rate. The PS280 or PS283 then crosses over to constant-voltage mode.


Figure 3: Constant Voltage/Constant Current Crossover

## Setting the Current Limit

Before you begin using the PS280 or PS283 to power a device, you should set its current limit lower than the maximum safe current for the device to be powered.

**CAUTION.** In order to avoid damaging your device with a current overload, set the current limit on the PS280 or PS283 before you connect it to your device.

1. Determine the maximum safe current for the device to be powered.
2. With the test lead, temporarily short the positive and the negative output terminals of the power supply together.

3. Rotate the **VOLTAGE** knob away from zero sufficiently to light the **C.C.** indicator.
4. Set the meter selection switch to **AMPS** so that the LED display shows the current.
5. Adjust the **CURRENT** knob for the desired current limit.
6. Read the value shown on the LED display. This is your preset current limit. Do not increase the current control setting.
7. Remove the short between the positive and negative output terminals.

You are now ready to power your device.

## Test Modes

The two variable power supplies on your PS280 or PS283 can be operated independently of each other, or the slave supply can track the master supply. Below are instructions for operating the instrument in independent modes, followed by instructions for operating the instrument in series or parallel tracking modes.

### Independent Modes

In independent mode, any one output of each power supply can be connected to any one terminal of another supply or to ground. The variable supplies are independently controlled by the front panel **VOLTAGE** and **CURRENT** control knobs.

There are three independent modes in which you can operate the PS280 or PS283: floating, ground-referenced, and stacked.


In floating mode, the power supply is not referenced with respect to ground.

In ground-referenced mode, one of the output terminals is grounded, providing a fixed reference point for your measurement.

In stacked mode, you connect the negative output terminal of one variable power supply to the positive output terminal of the other. The stacked configuration allows you to test a circuit requiring between 30 and 60 V. A stacked configuration can be either floating or ground-referenced.

**Floating.** In the independently floating mode, each variable power supply provides from 0 to 30 V at 0 to 2 A (0 to 1 A for the PS283).

Figure 4 shows each of the three power supplies connected to a separate load.

**Figure 4: Independent Floating Application**

The tracking switches are disengaged for independent operation. The left voltage and current control knobs control the outputs for the slave variable power supply, and the right knobs do the same for the master power supply. All outputs are electrically independent.

To test a circuit in the independently floating mode, follow these steps:

1. Press the **POWER** button to apply power to the PS280 or PS283.
2. Rotate the **VOLTAGE** knob to zero.
3. Determine the polarity of your device.


4. Plug one of the test leads into the positive output terminal.
  5. Plug the other test lead into the negative output terminal.
  6. Press **POWER** to turn off the PS280 or PS283.
  7. Clip the positive test lead to the positive pole of your device.
  8. Clip the negative test lead to the negative pole of your device.
  9. Press **POWER** to turn on the PS280 or PS283.
- 10.** Push the **AMPS/VOLTS** selection switch so that the LED display shows either voltage or current, as you want.
- 11.** Rotate the **VOLTAGE** knob as desired.
- 12.** If you are using a preset current limit (see page 10), do not touch the **CURRENT** knobs. Otherwise, rotate the **CURRENT** knob as desired.

**Ground-Referenced.** In the independently ground-referenced mode, each variable power supply provides from 0 to 30 V referenced with respect to ground at 0 to 2 A (0 to 1 A for the PS283). Any one of a pair of output terminals, either the positive or the negative, can be connected to ground. The FIXED 5 V power supply can also be ground-referenced.

**WARNING.** *In order to avoid grounding the power line, which can cause electrical shock, explosion, or fire, isolate the device being powered from the line voltage power source when using any ground-referenced output configuration from the PS280 or PS283.*

Figure 5 shows an example of a circuit with the FIXED 5 V terminal referenced to ground and both the master and slave variable power supplies referenced to -5 V.

**Figure 5: Independent Common Ground-Referenced Application**

In this configuration, each of the variable power supplies can be varied from  $-5\text{ V}$  to  $+25\text{ V}$  ( $+30\text{ V}$  overall). The GND post becomes the relative negative terminal for both variable outputs. Because the variable power supplies are referenced to  $-5\text{ V}$ , the LED display, when set to display volts, shows a value that is five volts lower than the actual output.

For example, the LED display indicates:

- $0\text{ V}$  when the output is  $-5\text{ V}$
- $5\text{ V}$  when the output is  $0\text{ V}$
- $30\text{ V}$  when the output is  $25\text{ V}$

Negative  $5\text{ V}$  is available between GND and the negative terminal of the FIXED  $5\text{ V}$  power supply.

To test a circuit in the independently ground-referenced mode, follow these steps:

1. Turn the **POWER** off to the PS280 or PS283.
2. Connect the outputs as shown in Figure 5.
3. Set both variable supply **VOLTAGE** controls to the minimum setting.

4. Set both variable supply **CURRENT** controls to midrange.
5. Set the **AMPS/VOLTS** switches for both power supplies to display volts.
6. Turn on the **POWER** to the PS280 or PS283. The display should read 0 V for both variable power supplies. An external meter connected across the load or load terminals should read -5 V.
7. Turn the **POWER** off to the PS280 or PS283 again.
8. Connect the device or devices to be tested.
9. Turn on the **POWER** to the PS280 or PS283 again. Adjust the voltages as needed.

Figure 6 shows the PS280 or PS283 connected to produce separate outputs of +5 V from the FIXED power supply, 0 to +30 V from the slave variable power supply, and 0 to -30 V from the master variable power supply. In this configuration, the red output terminal of the master variable power supply is the negative reference terminal because it is directly connected to the ground terminal.



Figure 6: Independent Ground-Referenced Split Application

To test a circuit in an independent ground-referenced split application mode, follow these steps:


1. Turn the **POWER** off to the PS280 or PS283.
2. Connect the outputs as shown in Figure 6.
3. Set both variable supply **VOLTAGE** controls to the minimum setting.
4. Set both variable supply **CURRENT** controls to midrange.
5. Turn on the **POWER** to the PS280 or PS283.
6. Set the desired voltages for both variable power supplies.
7. Turn the **POWER** off to the PS280 or PS283 again.
8. Connect the device or devices to be tested.
9. Turn on the **POWER** to the PS280 or PS283 again. If necessary, readjust the voltages.

Figure 7 shows the configuration for three ground-referenced negative power supplies.

**Figure 7: Three Ground-Referenced Negative Power Supplies**


Figure 8 shows the configuration for three ground-referenced positive power supplies.

**Figure 8: Three Ground-Referenced Positive Power Supplies**

**Stacked.** In the independently stacked mode, the variable power supplies are connected and provide from 0 to 60 V at 0 to 2 A (0 to 1 A for the PS283).

Figure 9 on page 18 shows the PS280 or PS283 connected in a stacked manner to produce a variable output of 0 to +60 V ground-referenced. The FIXED power supply produces -5 V, ground-referenced. In this configuration, the red output terminal of the master variable power supply is the negative reference terminal because it is directly connected to the ground terminal.

**Figure 9: Independent Positive Stacked Application**

To test a circuit in the independently stacked mode, follow these steps:

1. Turn the **POWER** off to the PS280 or PS283.
2. Connect the outputs as shown in Figure 9.
3. Set both variable supply **VOLTAGE** controls to the minimum setting.
4. Set both variable supply **CURRENT** controls to midrange.
5. Turn on the **POWER** to the PS280 or PS283.
6. Set the desired voltage. Observe the LED display; the total output is the sum of both voltage readings.
7. Turn the **POWER** off to the PS280 or PS283 again.
8. Connect the device or devices to be tested.
9. Turn on the **POWER** to the PS280 or PS283 again. If necessary, readjust the voltages.

Figure 10 shows a stacked application in which you have a 0 to -60 V output from the variable power supplies and a +5 V output from the FIXED power supply.



Figure 10: Independent Negative Stacked Application

## Tracking Modes

There are two tracking modes in which you can operate the PS280 or PS283: series and parallel.

**Series.** In series mode, the positive output terminal of the master variable power supply is internally connected to the negative output terminal of the slave power supply. This connection allows the PS280 or PS283 to produce 0 to 60 V at 0 to 2 A (0 to 1 A for the PS283).

When you place the PS280 or PS283 in series mode, the output terminals are hooked together internally as shown in Figure 11.

**Figure 11: Series Tracking Inside the PS280 or PS283**

The voltage knob for the master variable power supply controls the voltage for both variable power supplies. Using the master voltage control, the maximum slave supply voltage is automatically set to the same value as the master supply.

To test a circuit in the series tracking mode, follow these steps:

1. Turn the **POWER** off to the PS280 or PS283.
2. Connect the outputs as shown in Figure 12.


Figure 12: Series Tracking Application

3. Set the PS280 or PS283 to series tracking mode by pressing the left **TRACKING** button. Make sure that the right **TRACKING** button is released (out).
4. Set the master **AMPS/VOLTS** switch to the voltage metering position. Set the slave **AMPS/VOLTS** switch to the current metering position. This allows you to simultaneously monitor both current and voltage.

**NOTE.** In series tracking mode, the output voltage is double the value displayed on the voltage metering LED display, because both supplies are producing the same voltage.

5. Set the slave **CURRENT** knob fully clockwise.
6. Set the current limit using the master **CURRENT** knob. (Refer to *Setting the Current Limit* on page 10.)

---

**NOTE.** In series tracking mode, the current flowing through the two supplies must be equal. Therefore, the maximum current limit is the lower of the values set by the two current control knobs.

---

7. Turn on the **POWER** to the PS280 or PS283.
8. Adjust the output voltage to the desired level using the master **VOLTAGE** knob.
9. Turn the **POWER** off to the PS280 or PS283 again.
10. Connect the device or devices to be tested.
11. Turn on the **POWER** to the PS280 or PS283 again. Readjust the voltages if necessary.

---


**NOTE.** The 5 V FIXED supply can be independently grounded or allowed to float.

---

**Parallel.** In parallel tracking mode, the positive output terminals of both variable power supplies are internally connected, and the negative output terminals of both variable power supplies are internally connected. These connections allow the PS280 or PS283 to produce 0 to 30 V at 0 to 4 A (0 to 2 A for the PS283).

When you place the PS280 or PS283 in parallel mode, the output terminals are hooked together internally as shown in Figure 13.

The master power supply's **VOLTAGE** and the **CURRENT** knobs control the voltage and current for both variable power supplies.

**Figure 13: Parallel Tracking Inside the PS280 or PS283**

To test a circuit in the parallel tracking mode, follow these steps:


1. Turn the **POWER** off to the PS280 or PS283.
2. Connect the outputs as shown in Figure 14 on page 24.
3. Set the PS280 or PS283 to parallel tracking mode by pressing both tracking buttons.
4. Set the master **AMPS/VOLTS** switch to the voltage metering position, and set the slave **AMPS/VOLTS** switch to the current metering position. This allows you to simultaneously monitor both current and voltage.

---

**NOTE.** In parallel tracking mode, the output current is double the value displayed on the current metering LED display, because both supplies are producing the same amount of current.

---

5. Turn on the **POWER** to the PS280 or PS283.
6. Set the current limit using the master **CURRENT** knob. (Refer to *Setting the Current Limit* on page 10.)

**Figure 14: Parallel Tracking Application**

7. Adjust the output voltage to the desired level using the master **VOLTAGE** knob.
8. Turn the **POWER** off to the PS280 or PS283 again.
9. Connect the positive polarity of the device being powered to the positive master terminal.
10. Connect the negative polarity of the device being powered to the negative master terminal.

**CAUTION.** *To prevent damage to the PS280 or PS283, do not attempt to obtain output simultaneously from both variable power supplies while in parallel tracking mode.*

**NOTE.** *The 5 V FIXED supply can be independently grounded or allowed to float.*

# Appendix A: Specifications

**Table 1: Physical Characteristics**

Dimension	Measurement
Width	255 mm (10.0 in)
Height	145 mm (5.7 in)
Depth	335 mm (13.2 in)
Weight	11.5 kg (25.4 lb) PS280 9.0 kg (19.9 lb) PS283

**Table 2: Environmental Characteristics**

Characteristic	Temperature	Relative Humidity
Storage	-10°C to +70°C	70%
Operating	0°C to 40°C	80%

**Table 3: Operational Characteristics**

Characteristic	Measurement
Outputs	Two 0 to 30 VDC, one 5 VDC
Voltage (5 V)	5.0 ±0.25 VDC at 3.0 A maximum foldback current limited
Voltage (0–30 V)	0–30 constant VDC at 2.0 A constant, maximum (PS280) or 1.0 A constant, maximum (PS283)
Line Regulation (5 V)	≤5 mV
Line Regulation (CV)	≤0.01% +3 mV PS280 ≤0.01% +5 mV PS283

**Table 3: Operational Characteristics (Cont.)**

<b>Characteristic</b>	<b>Measurement</b>
Line Regulation (CC)	$\leq 0.2\% + 3 \text{ mA}$
Load Regulation (5 V)	$\leq 0.2\%$
Load Regulation (CV)	$\leq 0.01\% + 3 \text{ mV}$ (rating current $\leq 3 \text{ A}$ ) $\leq 0.01\% + 5 \text{ mV}$ (rating current $> 3 \text{ A}$ ) $\leq 300 \text{ mV}$ (0–60 V single series tracking supply)
Load Regulation (CC)	$\leq 0.2\% + 3 \text{ mA}$
Ripple/Noise (5 V)	$\leq 2 \text{ mV rms}$
Ripple/Noise (CV)	$\leq 1 \text{ mV rms}$ , 5 Hz–1 MHz
Ripple (CC)	$\leq 3 \text{ mA rms}$
Temperature Coefficient (CV)	$\leq 300 \text{ ppm}/^\circ\text{C}$
Recovery Time (CV)	$\leq 100 \mu\text{s}$ (time to recover after a 50% load change with 0.5 A minimum)
Tracking Error (Slave)	$\leq 0.5\% + 10 \text{ mV}$ of the master supply
Indicator	Two 3 1/2 digit 0.5 in LED panel display meter
Meter Indicators	0–30 VDC $\pm(0.5\% \text{ of reading} + 2 \text{ digits})$ 0–2 A $\pm(0.5\% \text{ of reading} + 2 \text{ digits})$
Insulation (Chassis-to-Terminal)	$\geq 20 \text{ M}\Omega$ at DC 500 V
Insulation (Chassis-to-AC Cord))	$\geq 30 \text{ M}\Omega$ at DC 500 V

**Table 4: Electrical Characteristics**

<b>Characteristic</b>	<b>Measurement</b>
Line voltage	90 to 110 108 to 132 198 to 242 216 to 250, all VAC at 50–60 Hz
Power consumption	386 VA, 300 W maximum (PS280) 265 VA, 200 W maximum (PS283)

**Table 5: Certifications and Compliances**

EC Declaration of Conformity – EMC	<p>Meets intent of Directive 89/336/EEC for Electromagnetic Compatibility. Compliance was demonstrated to the following specifications as listed in the Official Journal of the European Communities:</p> <table> <tbody> <tr> <td>EN 55011</td><td>Class B Radiated and Conducted Emissions</td></tr> <tr> <td>EN 50081-1 Emissions:</td><td></td></tr> <tr> <td>    EN 60555-2</td><td>AC Power Line Harmonic Emissions</td></tr> <tr> <td>EN 50082-1 Immunity:</td><td></td></tr> <tr> <td>    IEC 801-2</td><td>Electrostatic Discharge Immunity</td></tr> <tr> <td>    IEC 801-3</td><td>RF Electromagnetic Field Immunity</td></tr> <tr> <td>    IEC 801-4</td><td>Electrical Fast Transient/Burst Immunity</td></tr> <tr> <td>    IEC 801-5</td><td>Power Line Surge Immunity</td></tr> </tbody> </table>	EN 55011	Class B Radiated and Conducted Emissions	EN 50081-1 Emissions:		EN 60555-2	AC Power Line Harmonic Emissions	EN 50082-1 Immunity:		IEC 801-2	Electrostatic Discharge Immunity	IEC 801-3	RF Electromagnetic Field Immunity	IEC 801-4	Electrical Fast Transient/Burst Immunity	IEC 801-5	Power Line Surge Immunity
EN 55011	Class B Radiated and Conducted Emissions																
EN 50081-1 Emissions:																	
EN 60555-2	AC Power Line Harmonic Emissions																
EN 50082-1 Immunity:																	
IEC 801-2	Electrostatic Discharge Immunity																
IEC 801-3	RF Electromagnetic Field Immunity																
IEC 801-4	Electrical Fast Transient/Burst Immunity																
IEC 801-5	Power Line Surge Immunity																
EC Declaration of Conformity – Low Voltage	<p>Compliance was demonstrated to the following specification as listed in the Official Journal of the European Communities:</p> <p>Low Voltage Directive 73/23/EEC, amended by 93/68/EEC.</p> <table> <tbody> <tr> <td>HD401 S1</td> <td>Safety Requirements for Electronic Measuring Apparatus.</td> </tr> </tbody> </table>	HD401 S1	Safety Requirements for Electronic Measuring Apparatus.														
HD401 S1	Safety Requirements for Electronic Measuring Apparatus.																

## **Appendix B: Maintenance**

This appendix provides information for the basic maintenance of the PS280 or PS283 Laboratory DC Power Supply.

### **Cleaning**

To clean the Laboratory DC Power Supply, use a soft cloth dampened in a solution of mild detergent and water. Do not spray cleaner directly onto the instrument, since it may leak into the cabinet and cause damage.

Do not use chemicals containing benzine, benzene, toluene, xylene, acetone, or similar solvents.

Do not use abrasive cleaners on any portion of the power supply.

### **Preparing for Shipment**

If the original packaging is unfit for use or not available, use the following packaging guidelines:

1. Use a corrugated cardboard shipping carton having inside dimensions at least three inches greater than the instrument dimensions.
2. Put the instrument into a plastic bag or wrap to protect it from dampness and loose packing material.
3. Place the instrument into the box and firmly stabilize it with packing material.
4. Seal the carton with shipping tape.

## Troubleshooting

Electronic maintenance on the power supply must be performed by a trained technician. However, an operator can perform some basic and routine maintenance. Perform the following steps to isolate the fault:

1. The power switch is on. The instrument is plugged in. Neither the C.C. nor the C.V. indicator is lighted.

Check the output terminals with a voltmeter.

- a. Set the voltage control of the voltmeter to midrange.
- b. Ensure that the range and polarity settings are correct.
- c. Place the voltmeter jacks in the PS280 or PS283 output terminals.
- d. Determine if the terminals are producing any output.

Are the outputs working?

Yes                    Go to step 2.

No                    Go to step 3.

2. Refer to a service technician.

***WARNING.*** To prevent electrical shock, unplug the power cord and disconnect the test cables from any power source before checking or replacing the fuse.

---

3. Check the fuse with a multimeter.

- a. Set the multimeter to the low ohms range.
- b. Apply the multimeter probes across the fuse.
- c. Determine if a continuous circuit exists.

Is the fuse okay?

Yes                    Go to step 5.

No                    Go to step 4.

4. Replace the fuse.

5. Verify that the line settings on the rear panel match the line voltage. Do they?  
Yes                  Go to step 7.  
No                  Go to step 6.
6. Reset the line settings. Refer to *Preparing the Power Supply for Use* on page 2.
7. Check the power cord.

***WARNING.*** To prevent personal injury, be sure the power cord is disconnected at both ends before you check it.

- Is the power cord frayed or broken?
- Yes                  Go to step 8.  
No                  Go to step 2.
  8. Replace the power cord.
  9. The power switch is on. A variable output power supply is connected to a circuit. The C.C. or C.V. indicator is on. Neither variable output power supply is producing any electrical output.  
Disconnect the instrument from the circuit. Check the output terminals with a voltmeter. Are the outputs working?  
Yes                  Go to step 10.  
No                  Go to step 2.
  10. Check the circuit you have been testing for a short or low resistance.

## Appendix C: Replaceable Parts

Replaceable parts may be ordered directly from your authorized Tektronix dealer.

### Standard Accessories

The following items are shipped with the Laboratory DC Power Supply:

**Table 6: Standard Accessories**

Accessory	Tektronix Part Number
Fuse, 5 x 20 mm, 4 A, 250 V, SB (PS280: 90 – 132 V operation)	159-0297-00
Fuse, 5 x 20 mm, 2.5 A, 250 V, SB (PS283: 90 – 132 V operation)	159-0226-00
Test Leads	196-3201-00
User Manual	070-8355-XX
115 V Power Cord	Refer to Table 8

### Optional Accessories






The following items are available as optional accessories:

**Table 7: Optional Accessories**

Accessory	Tektronix Part Number
Fuse, 5 x 20 mm, 2 A, 250 V, SB (PS280: 198 – 250 V operation)	159-0107-00
Fuse, 5 x 20 mm, 1.25 A, 250 V, SB (PS283: 198 – 250 V operation)	159-0247-00
230 V Power Cords	Refer to Table 8

The following power cords are available:

**Table 8: Accessory Power Cords**

Plug Configuration	Normal Usage	Tektronix Part Number
	North America 115 V	161-0104-00
	Europe 230 V	161-0104-06
	United Kingdom 230 V	161-0104-07
	Australia 230 V	161-0104-05
	North America 230 V	161-0104-08
	Switzerland 230 V	161-0167-00